Introduction

Congratulations on your purchase of LEDSONIC+ Blue LED. We welcome you to experience the phototherapy effects of LEDSONIC+ complete system! Please check the attachment carefully upon purchase and read this manual before use, ensuring that you completely understand the correct use of this product. LEDSONIC+ is a hand held device, allowing professional acne treatment in the convenience of your home. The product's blue light treatment has a therapeutic effect on acne while the red light treats the inflammation caused by acne.

How does LEDSONIC+ LED works?

For a long time, blue light demonstrated its antibacterial properties; specific wavelengths of blue light (415nm) penetrate the skin and help destroy acne-causing bacteria. It reduces inflammatory acne lesions, bringing an effective and safe method in the treatment of acne.

According to the nature of skin, the severity of acne can be divided into 4 grades: grade 1 (mild): only acne; grade 2 (moderate): acne and inflammatory papules; grade 3 (severe): in addition to acne, the skin has inflammatory papules and pustules; grade 4 (very severe): in addition to acne, inflammatory papules and pustules, nodules, cysts or scars appear. Acne classification is an important basis of acne treatment, and curative effect. LEDSONIC+ therapeutically assists in mild-to-moderate acne cases. As it can be used on the face it is especially effective for improving its appearance. The treatment head can only be used on the forehead, face and under eyes. The device's head must be kept a small distance away from eyes (1 cm.) Do not point directly into the eyes, and direct contact with the neck is not allowed.

If facial skin is inflamed with redness, etc., please avoid direct contact with the treated skin.

During treatment process, the product can be used twice a day: 3 minutes of blue light. To achieve a greater results, we recommend you to use with LEDSONIC+ red light (not included in this package) for anti-inflammatory and healing. It reduces inflammatory acne lesions, bringing an effective and safe method in the treatment of acne. After used 3 minutes of blue light, we recommend to follow by 3 minutes of red light, with the device in close contact on acne cysts. With continuous usage, you may find a decrease in acne inflammation, acne papules, and pain, for an acne free, perfect complexion.

Product specific use:

The product's intended use: The treatment of facial acne and elimination of inflammation caused by facial acne. The product's blue-ray treatment head mainly has a therapeutic effect on facial acne, with greater results when used together with red light. Cosmetically, it diminishes wrinkles, minimizes pores, improves skin texture and color, and improves dark under eyes circles. In addition it reduces scars, freckles, local pigmentation, and skin dehydration.

Disclaimer:

- When used according to instructions, side effects and complications associated with use of LEDSONIC+ are rare. However, every energy based system may involve some degree of risk. Therefore, it is important that you understand and accept the risks and complications that can occur. LEDSONIC+ is a gentle home use device and most users only sense moderate heat during use with no further discomfort.
- Do not use: If you are taking photosensitive drugs; Have a thyroid disorder and medication; If the acne has reached level 3; If you have been injected with steroids; Treatment of parts with open wounds; Face treatments with liniments, salves or ointments that contain heat - producing ingredients.
- Do not use if you are pregnant, breast feeding or nursing.
- Do not use if you have a pacemaker or internal defibrillator.
- Do not use if you have superficial metal or other implants in the desired treatment area.
- Do not use on dry or fragile skin caused by the use of chemical peels, glycolic peels, Alpha Hydroxy Acids (AHAs), Beta Hydroxy Acids (BHAs), Retin-A, Accutane and topical isotretinoin and azelaic acid.
- Do not use if you have varicose veins in the treatment area.
- Do not use if you have eczema, psoriasis, lesions, open wounds or active infections, such as cold sores in the area to be treated. Wait for the affected area to heal before using LEDSONIC+.
- Do not use if you have an impaired immune system due to immune suppressive diseases such as AIDS and HIV, or use of immunosuppressive medications.
- Do not use if your skin is sensitive to light, easily developed a rash or an allergic reaction.
- Do not use if you have had surgery in the areas to be treated.
- Do not use if you have diabetes, or other systemic or metabolic diseases or

congestive heart disease.

- Do not use if you have epilepsy.
- Do not use if you have any other condition which in your physician's opinion would make it unsafe for you to be treated.
- Do not service /maintain the device whilst in use.
- Keep away from children.
- If you have any problems with this device, such as setting up, maintaining or usage, please contact your local provider or contact us at www.qlarite.com. Don't open or repair the device by yourself.
- Please report to your local provider if any unexpected operation or event occurs.
- When charging/using, do not touch charging connector and the patient simultaneously. Charger/charger base is specified as a part of the EQUIPMENT.

Package contents:

- Blue Light Head (wavelength: 415nm)
- Travel Pouch

Guarantee

We warrant this product against any defects that are due to faulty material or workmanship for a 2 year period from the original date of consumer purchase. If the product should become defective within the warranty period, we will repair any such defect without charge provided there is proof of purchase. This does not mean an extension of warranty period. The warranty shall apply in all countries in which our product was sold via an authorized dealer.

This warranty does not cover the damage is due to by accident or misuse, abuse, alteration to the product or use inconsistent with the technical and /or safety instructions required.

This warranty shall not apply if the product has been dismantled or repaired by a person not authorized by us.

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Copyright © 2017 Qlarité SAS. All Rights Reserved.
Qlarité SAS