LEDSONIC+

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ENG

LEDSONIC+

All in One Skin Cleansing and Anti aging system:

- Deep skin cleansing
- Increases skin radiance, density
- Evens skin tone
- Reduces wrinkles and fine lines
- Improves acne-prone skin

User Guide

Important Safety Warnings!

Do not stare into the light!

When using the LEDs attachments, please do not stare directly into the light as this will cause discomfort or injury.

When using the LEDs attachments, do not use on your eyes or eyelids!

When using the LEDs attachments, do not use over the thyroid gland (neck area).

When using the ultrasonic attachments, keep away from eyes by at least 1cm away at all times!

Caution: Protect your eyes from blue or red light.

Children should be supervised at all times to ensure they do not play with the device.

Please read all warnings and safety precautions before use, and strictly follow the guidelines when using LEDSONIC+.

Keep away from the children. Close supervision is necessary when this product is used by, on, or near children or physically challenged persons.

Not servicing/maintenance while the me equipment is in use.

Table of Contents

1. Introduction
2. Benefits
3. Before You Start
 - 3.1 General Description
 - 3.2 Package contents
4. Important Safety Precautions
 - 4.1 Warning
 - 4.2 Precautions
 - 4.2.1 Avoid to use LEDSONIC+ if...
 - 4.2.2 Areas not to treat
 - 4.3 Possible side effects
5. Preparing for use
 - 5.1 Charging
 - 5.2 Preparing treatment areas
 - 5.3 Patch test
6. Using LEDSONIC+
 - 6.1 LEDSONIC+ setup
 - 6.1.1 How to change attachment heads
 - 6.1.2 Operating the system
 - A) Ultrasonic therapy mode
 - B) Light therapy mode
 - C) Using sonic cleansing brush
 - 6.2 Treatment tips
 - A) Ultrasonic therapy mode
 - B) Light therapy mode
 - C) Sonic cleansing brush mode
7. Cleaning, maintenance and storage
 - 7.1 Cleaning
 - 7.2 Maintenance and storage
 - 7.3 Protecting the environment
8. Troubleshooting
9. Specification
10. Labeling
11. Service and Warranty
12. Copyright Reserved

ENG

ENG

1. Introduction

Congratulations on your purchase of LEDSONIC+! Experience the 4 in 1 therapy effects of LEDSONIC+!

How does Ultrasonic works?

The application of ultrasonic frequency is commonly used by beauty and clinic professionals but now it is available to home users with LEDSONIC+. It offers safe and effective facial toning treatments in the comfort of your own home. It plumps up and smoothes fine lines and wrinkles, tightens skin, shrinks enlarged pores, relieves puffy eyes, improves dark eye circles, stimulates cellular metabolism, improves and reduces scars, freckles, age spots and redness. The process is safe, painless and gentle.

What is the difference between 1MHz and 5MHz?

Generally, lower ultrasonic frequencies will penetrate deep below the skin. Thus, the 1MHz mode is designed for slightly thicker facial tissues and muscles (i.e. cheeks, chinetc.). The 5MHz ultrasonic mode is appropriate for treatment around the most delicate and thinner areas (i.e. under eyes, temples, forehead and sides of the nose)

How does the ultrasonic treatment offer improvement of your skin?

When ultrasounds are absorbed by muscles, tendons and bones, the harder surfaces like bone absorb and reflect ultrasound energy. This in turn dramatically increases oxygen supply to the cells. This antibacterial cleansing action also helps to improve acne and blemishes. The skin is left feeling instantly energized, firm, tighter, refreshed and noticeably softer – even after just one application. (See Fig.1) The intensive rejuvenating cyclic soundwaves penetrate below the skin, reaching the dermis. Generally, 1MHz frequency mode will penetrate deeper below the skin. It is designed for the face (except around the eye area) where tissues and muscles are much thinner. The 5MHz frequency mode is appropriate for treatment around the eyes, which is the most delicate area. When used with skin care products, ultrasounds can help the ingredients penetrate deep into the skin. It can maximize the effects of your skin care products.


Fig. 1

Stimulated by light and ultrasound, the skin becomes firmer and more elastic

How does light therapy by LEDSONIC+ work?

Visible light therapy, conceived at the end of the 19th century by Niels Fynsen, a Danish doctor, won him the Nobel Prize of Medicine in 1903. Almost a century later, thanks to space research, the amazing biological properties of light and their stimulation of the cells' respiratory system were finally brought to light. These properties manifest themselves especially with red and infrared, as well as blue. They activate repair and healing of scars, fight inflammation in tissues. Numerous publications illustrate these properties, and some of the most important studies come from research conducted by NASA. Today, LED technology can manufacture powerful light sources in the red and blue range. Red LED therapy by LEDSONIC+ also has an important role to play in skincare, notably with problem skins and acne. This indication is treated with Blue light, and was approved by the FDA in August 2002.

How does Red light therapy by LEDSONIC+ benefit skin?

LEDSONIC+ uses a super-luminous LED array to deliver concentrated beams of red light (660nm). The phototherapy enhances blood circulation and stimulates fibroblast collagen production to strengthen dull and wrinkled skin. This cellular and metabolic effect helps fight skin aging, improves skin tone and increases tissue oxygenation. Red light therapy is safe for all skin types and ages, there is no down time or

recovery time.

It is one of the most effective and non-invasive ways to help rejuvenate and maintain youthful skin.

How does the sonic cleansing brush by LEDSONIC+ give you a smoother and younger looking skin?

LEDSONIC+ cleansing brush provides a much deeper cleanse than washing by hand to help gently exfoliate dead skin cells and remove dirt and oil. The silicone extra hygienic and ultra-gentle brush delivers over 8000 times sonic pulsations to deep cleanse effectively, enhancing the absorption of your favorite skincare products.

2. Benefits

LEDSONIC+ provides you with a total skin care solution:

Mode	Treatment
Ultrasonic - 1 MHz (2 Level)	For cheeks, chin, the lower jaw (See Fig.2 "U1")
5 MHz (2 Level)	For under and around eyes, forehead, temple, neck (See Fig.2 "U2")
Red light	Skin rejuvenation, skin tightening and skin care penetration (See Fig.2)
Sonic cleansing brush	Skin deep cleansing (See Fig.2)

Fig. 2

Advanced System formulated by leading dermatologists

Our PT3® patented technology, clinically proven in dermatology, restores radiance, reverses aging by improving the density of the skin and dramatically reduces wrinkles and fine lines.

Unique pattern design to maximize effectiveness

A unique pattern design on the surface of the ultrasonic attachment holds effectively the ingredients of your skin care products during contact with the skin. (See Fig.3)


Fig. 3

High intensity energy

High power density of 50mW/cm² illuminated surface produces 9 joules/cm² dose in 180 seconds.

Simple and easy to use

All in one anti aging and skin care device, easy to operate.

Gentle treatment

During the process you will only feel a little heat. It does not contain any harmful ultraviolet (UV) light which can damage skin and treatment does not result in dryness or discomfort.

Extra hygienic

LEDSONIC+ sonic cleansing brush is made of quick drying and impermeable silicone that resists bacteria buildup and doesn't require replacement brush heads.

ENG
ENG

3. Before You Start

3.1 General Description (Fig. 4)

1. Main Body
2. Attachment connector base
3. Touch key for "Power On"
4. Touch key for "Select" : Ultrasonic mode or light on
5. Ultrasonic mode indicator
6. Charge contact point
7. Charger Base
8. 1MHz/ 5MHz Ultrasonic attachment
9. Red LEDs attachment
10. Sonic cleansing brush attachment
11. Power adaptor
12. Silicon Protector


Fig. 4

Notes:

- Applied parts are included the ultrasonic head and red LED and sonic cleansing brush heads.
- The AC/DC adaptor is used to isolate the equipment from the supply mains.
- Charger/charger base is specified as part of the device.

3.2 Package contents

- Applicator
- Attachment for Ultrasonic/ Red LEDs/ Sonic cleansing brush
- Silicon Protector
- Charger Base
- Adaptor
- Instructions manual
- Soft Cloth

4. Important Safety Precautions

Please read this instruction manual carefully before you use the appliance and save it for future reference. If you are unsure if LEDSONIC+ is safe for you to use, please consult your doctor or dermatologist!

4.1 WARNING:

As with most electrical appliances, electrical parts are electrically live even when the switch is off. To reduce the risk of injury or death by electric shock:

- Do not reach for the unit if it has fallen into water or other liquid. Unplug unit immediately. Do not store charger where it can fall or be pulled into a sink or bathtub.
- Do not place or drop into water or other liquids.
- The adaptor and charger base are not waterproof, keep away from water.
- Keep the appliance and the adaptor dry at all times.
- Do not attempt to open or repair the appliance by yourself.
- If the appliance is broken, do not touch any inner part to avoid electric shock.
- Never insert any objects into the appliance.
- Do not use the appliance if it is not working properly or if it appears damaged.
- Do not use near the shortwave or microwave equipment.
- Make sure the adaptor does not get wet.
- Keep away from flammable aerosol products when used or where oxygen is administered.
- Only use parts supplied with the appliance.
- Only charge the appliance with the adaptor supplied.
- Do not use the appliance or the adaptor if it is damaged.
- Do not charge when exposed to direct sunlight or heat.
- Keep this product out of reach of children. The use of this appliance by persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge can cause hazards, unless under supervision or instruction concerning the use of appliance by a person responsible for their safety.
- The adaptor contains a transformer. Do not cut the adaptor to replace with another plug as this will cause a hazardous situation.
- Turn off the appliance if it becomes unusually warm to touch.
- Switch off the appliance before cleaning.
- Stop immediately if you experience any adverse effects or develop prolonged redness. You may be experiencing a photosensitive reaction.
- LEDSONIC+ is designed for home use (as opposed to exhaustive commercial

ENG

ENG

or professional use). Avoid such exhaustive use.

- LEDSONIC+ is only designed for personal use, **DO NOT** share with others.
- The patient is an intended operator.
- Do not use any accessories, detachable parts and materials which are not described in the instructions for use.
- Do not replace the battery by yourself as it can result in an unacceptable risk.
- No modification of this equipment is allowed.
- Do not touch the output the AC/DC adaptor when using the equipment.
- Do not inhale or swallow the small parts.
- Do not interconnect this equipment with other equipment.

4.2 Precautions

Never use the appliance with the LED attachments around and on the eyes to avoid the risk of eye damage.

Please note that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

4.2.1 Avoid to use LEDSONIC+ if...

- Do not use if you are pregnant, breast feeding or nursing.
- If you have a pacemaker or internal defibrillator.
- If you have superficial metal or other implants in the desired treatment area.
- Do not use on dry or fragile skin caused by the use of chemical peels, glycolic peels, Alpha Hydroxy Acids (AHAs), Beta Hydroxy Acids (BHAs), Retin-A, Accutane and topical isotretinoin and azelaic acid.
- If you have a history of any type of cancer, or pre-malignant moles.
- If you have varicose veins in the treatment area.
- If you have pre-cancerous lesions or multiple atypical moles in the areas to be treated.
- If you have eczema, psoriasis, lesions, open wounds or active infections, such as cold sores in the area to be treated. Wait for the affected area to heal before using LEDSONIC+.
- If you have an impaired immune system due to immunosuppressive diseases such as AIDS and HIV, or use of immunosuppressive medications.
- If your skin is sensitive to light and easily develops a rash or an allergic reaction.

- If you have had surgery in the areas to be treated.
- If you epilepsy with light sensitivity.
- If you have diabetes, or other systemic or metabolic diseases or congestive heart disease.
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.
- If you have epilepsy.
- If you are on gold therapy.
- If you have received radiation therapy or chemotherapy treatments within the past 3 months.
- If the ultrasonic treatment head has cracks.

Note: This list is not exhaustive. If you are not sure whether you can use the appliance, please consult a doctor.

4.2.2 Areas not to treat

- Do not use directly on eyes. Keep away from eyes by at least 1cm.
- Do not use on the middle of neck (throat).
- Do not use on nipples, areola, or genitals.
- Do not use on an area of recent surgery, deep peel, laser resurfacing, scars, or skin that has been damaged with burns or scalds.
- Remove the applicator from the skin if it is unusually warm.

4.3 Possible side effects

When used according to instructions, side effects and complications associated with use of LEDSONIC+ are rare. However, every energy based system may involve some degree of risk. Therefore, it is important that you understand and accept the risks and complications that can occur. LEDSONIC+ is a gentle home use appliance, most users only sense moderate heat applied by the device during use with no further discomfort.

- Skin may become slightly irritated after using LEDSONIC+. This is a typical and normal phenomenon. If the irritation does not disappear within 3 days, please consult a physician or a dermatologist.
- Skin may become slightly red. This should disappear a few hours after treatment.

ENG

ENG

5. Preparing for use

5.1 Charging

- Fully charge the batteries before you use LEDSONIC+ for the first time and when the batteries are empty. Fully charging the batteries takes up to 4 hr. 30 mins
- Low battery alert: Charge LEDSONIC+ when hearing a continuous "D..D..D" sound. This indicates the battery is low on power.
- Fully charged batteries provide operation time at least approx. 1 hr. 20 mins. for ultrasonic mode and approx. 1hr. 40 mins. for LEDs light therapy mode.
- Fully charge the appliance every 3 to 4 months, even if you do not use LEDSONIC+ for a long time.

Directions for charging LEDSONIC+:

1. Switch off the appliance.
2. Insert the small plug from the adaptor into the back of the charger base and put the adaptor in the wall socket. (See Fig. 5)


Fig. 5

Notes: The LEDs attachments in the charger base are only for storage purpose.

3. Charging indicator light (See Fig. 6)
 - The "SELECT" key blinks orange whilst charging.
 - The "SELECT" key illuminates in continuous orange when fully charged.


Fig. 6

Notes: It is normal when the adaptor and the appliance feel warm during charging.

Warning: NEVER cover the appliance and adaptor during charging. When charging/using, do not touch charging connector and any person simultaneously. Once device is fully charged, please disconnect. Do not overcharge.

4. Remove the adaptor from the wall socket and pull the small plug out of the appliance when fully charged.

5.2 Preparing treatment areas

- Start by cleaning your face with a non-oil based cleanser.
- Use the provided soft cloth to clean the surface of the attachment.

5.3 Patch test

- Familiarize yourself with the features of the LEDSONIC+ device.
- Turn on the LEDSONIC+ with the ultrasonic mode and test on a small patch of skin and wait 24 hours.
- Light Test : Turn on the LEDSONIC+ with the Red LED attachment (and test again with optional Blue LED attachment, this attachment does not include in this package) and test on a small patch of skin with 3 minute treatment on the same area. Wait 24 hours.
- The patch test is to ensure there are no adverse reactions. If irritation occurs do not continue treatment and please consult a doctor or dermatologist.

6. Using LEDSONIC+

6.1 LEDSONIC+ setup

6.1.1 How to change attachment heads?

Press the buttons on the left and right sides and place the selected attachment onto the attachment connector base. (See Fig. 7)


Fig. 7

6.1.2 Operating the system

(Please check with section 3.1 "LEDSONIC+ device description")

A) Ultrasonic Therapy Mode

- Step 1: Ensure Ultrasonic attachment is correctly attached on the connector base.
- Step 2: Long touch "⏻" for power on.
- Step 3: Touch "SELECT" each time for selecting different frequency modes (1MHz/ 5MHz) and energy levels (Normal/ Strong).(See Fig. 8 for the sequences of the default setting when press "SELECT")


Fig. 8

Note: For details, please see section 6.2 "Treatment Tips"

- Step 4: To change the frequency mode or energy level, touch "SELECT" again. The light indicator (See Fig. 6) will blink in different status in each selection :

Frequency Mode	Level	Indicator Status	Operation time (auto power off)
1 MHz	Normal	"1MHz" Blinking slowly	15 mins
	Strong	"1MHz" Blinking quickly	10 mins
5 MHz	Normal	"5MHz" Blinking slowly	15 mins
	Strong	"5MHz" Blinking quickly	10 mins

- Step 5: Long touch “⏻” for power off when finished treatment.
- Step 6: Use a dry soft cloth to clean the surface of the attachment.
- Step 7: The appliance will switch off after 15 min.

Note: 10 min. is the recommended treatment time per session. The appliance will automatically power off after 15 min.

B) Light Therapy Mode

⚠ Caution: Protect your eyes from red light.

If you have to treat areas around the eyes with blue or red attachment, use the silicon protector provided with the device. Always keep away from eyes by at least 1cm.

- Step 1: Red LED attachment is correctly placed on the attachment connector base.
- Step 2: Long touch “⏻” for power on.
- Step 3: Touch “SELECT” for light on.
- Step 4: Long touch “⏻” for power off after completing treatment.
- Step 5: Use a dry soft cloth to clean the surface of the attachment.
- Step 6: LEDSONIC+ will “beep” for each minute during treatment time.
- Step 7: LEDSONIC+ will power off after 9 min.

Note: The appliance will automatically power off after 9 min. Please wait 8 min. for cool down before powering on again.

C) Using sonic cleansing brush

⚠ Caution: Never use cleaning products containing alcohol, petrol or acetone as they may irritate the skin.

- Step 1: Ensure sonic cleansing brush attachment is correctly attached on the connector base.
- Step 2: Remove all makeup, dampen face with water and apply your regular cleansing product onto the brush attachment.
- Step 3: Dampen face with water and apply cleansing lotion onto brush. Switch on device.
- Step 4: Long touch “⏻” for power on.
- Step 5: Touch “SELECT” to start
- Step 6: Start cleansing at your chin with gentle upward strokes along your jaw to your ear, moving gradually higher to cleanse the whole cheek.
- Step 7: Change to the other side of your face and forehead by working in circular motion. Lastly, cleanse both sides of the nose by gently gliding the brush head up and down.

- Step 8: Detach the sonic cleansing brush head to wash the brush surfaces with clean water and soap. Pat dry with a lint-free cloth or towel.

Note: 1~2 mins. is the recommended cleansing routine. The appliance will “beep” for each minute during cleansing mode.

6.2 Treatment Tips

A) Ultrasonic Therapy Mode

- Apply the skin care products (serum/lotion) on the ultrasonic attachment surface (See Fig. 9) and choose frequency according to treatment area (See Fig 2).

Fig. 9

Fig. 10

- Maximizing the effects of your skin care products – Ultrasonic waves are effective for product penetration with specific frequencies for different areas:
 - A) Use with products like cream, serums or lotion: Repeat each step a total of 5 times: glide LEDSONIC+ very slowly with a slight pressure. (See Fig. 10)
 - B) Use with products like facial mask: Repeat each step a total of 10 times: glide LEDSONIC+ very slowly with slight pressure. (See Fig. 10)

B) Light Therapy Mode

LEDSONIC+ uses a super-luminous LED array to deliver concentrated beams of light for phototherapy effect. This safe, painless, natural healing process converts light energy within skin cells, very similar to the way photosynthesis takes sunlight and turns it into food energy in plants. To maximize the effect, we highly recommend:

- For acne treatment: Use the Blue LED light therapy (Blue LED head is sold separately) for 3 minutes over affected areas twice daily. To maximize the results, follow with the red LED light therapy for 3 minutes over affected areas twice daily.
- Skin rejuvenation and firming: Use the red LED light therapy, holding the appliance for 3 minutes on different areas of the face (see section 4.2.2 “Areas not to treat”), then move to another area of the face. The appliance will “Beep” to remind you of every minute of treatment.

- We recommend to use the light therapy 3-4 times per week.
- If an area requires more rejuvenation, we suggested to extend the process slightly over 3 minutes in that area.

C) Sonic cleansing brush mode

Avoid using clay-based, silicone-based or grainy cleansers, exfoliators or scrubs as they may damage the soft silicone touch-points.

7. Cleaning, maintenance and storage

7.1 Cleaning

 **Caution:** Before cleaning your LEDSONIC+, make sure that the power switch is OFF and the power cord is disconnected from the charger base.

- Never clean LEDSONIC+ or any of its parts under the tap or in the dishwasher.
- Never use abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the appliance.
- Use the soft cloth included with your LEDSONIC+ appliance to clean the following parts:
 - The surface of the main body and charger base
 - The surface of the attachments
- Always clean the sonic cleansing brush after use with water and soap, pat dry with a lint-free cloth or towel.
- Avoid leaving your sonic cleansing brush in direct sunlight and never expose it to extreme heat or boiling water. The life span will greatly increase with proper maintenance and care.

7.2 Maintenance and storage

- Use the dry soft cloth, dip it in clear water to clean the surface of the ultrasonic head and LED light head after each application.
- Store LEDSONIC+ in a dry place at a temperature between -10°C and 45°C.
- If you want to replace a part or want to purchase an additional part, please go to your local provider or contact us at www.qlarite.com.

7.3 Protecting the environment

- At the end of its life, **DO NOT** discard LEDSONIC+ with normal household waste. Dispose of it at an official recycling collection point.

ENG

ENG

8. Troubleshooting

Refer to this troubleshooting guide if you experience any problems with LEDSONIC+, as this section addresses the most common problems you could encounter with this appliance. If you still have problems, please contact your local provider or contact us at www.qlarite.com. **Don't** open or repair the device by yourself.

Problem	Possible reason	Solution
LEDSONIC+ does not work	The rechargeable batteries are low in power The appliance is defective	Charge the appliance (see section 5.1 "Charging") Contact the nearest Qlarité dealer
LEDSONIC+ suddenly switches off	The rechargeable batteries are low in power	Charge the appliance
The appliance temperature rises during operation and recharge	This is normal	No action is required
Battery charging indicator light ("SELECT" key orange light) does not turn on when connecting the PLUG into the CHARGER BASE SOCKET	The POWER ADAPTOR is not plugged in the electrical outlet	Insert the PLUG into the CHARGER BASE SOCKET. Connect another wall socket to check if the electrical outlet is live If it is still does not charge, contact your nearest Qlarité dealer


9. Specification

Mode No.		QM-EF18
Software version		QM-EF18-SW1.0
Ultrasonic Frequencies		1MHz/ 5MHz
LEDs Light wavelength		660nm
LEDs Light Energy		Max. 50mW/ cm²
Lithium-ion battery		11.1V / 750mAh
Charging time		Approx. 4 hr. 30 mins
Battery lifetime		300 discharge cycles
Operating time	Ultrasonic Therapy	1 hr. 20 mins
	LEDs Light Therapy	1 hr. 40 mins
Adaptor spec.		Input: AC100-240V, 50-60Hz, 0.2A; Output: DC12V, 400mA Fuse(3N): AC250V; 1.0A, 3.6x10mm;Time-lag, low breaking capacity
Operating conditions		Temp. +15°C to +30°C
		Relative humidity: 25% to 75%
		Barometric pressure: 800-1060hPa
Storage and transport conditions		Temp. -10°C to +45°C
		Relative humidity: 10% to 90%
		Barometric pressure: 800-1060hPa

ENG

ENG

10. Labeling

Symbol	Description
	Warning/ Cautions
	Keep away from eyes
	CE Mark
	FCC Mark
	WEEE Waste Electrical and Electronic Equipment
	Refer to instruction manual
IPX5	Water projected by a nozzle (6.3 mm) against enclosure from any direction shall have no harmful effects. Test duration: 1 minute per square meter for at least 15 minutes Water volume: 12.5 litres per minute Pressure: 30 kPa at distance of 3 m (Only applied to Sonic cleansing brush head)
IPX4	Water splashing against the enclosure from any direction shall have no harmful effect, utilizing either: a) an oscillating fixture, or b) A spray nozzle with no shield. Test a) is conducted for 10 minutes. Test b) is conducted (without shield) for 5 minutes minimum. (Applied to LEDSONIC+ main device)
	Fragile; handle with care
	Keep the product in the dry place away from water and rain

Symbol	Description
	This way up
	Product package should be recycled

Safety Standards

- LVD: EN 60335-2-23:2003 + A1:2008 + A11:2010+A2:2015
EN 60335-1: 2012 + A11: 2014
EN 62233: 2008
- EMC: EN55011:2016, EN55014-2:2015, EN61000-3-2:2014,
EN61000-3-3:2013
- FCC: ANSI C63.4 & 47 CFR PART15, FCC OST/MP-5 & FCC PART 18
- Manufacturer System: ISO13485:2003 for Medical CE

ENG

ENG

11. Service and Warranty

We warrant this product against any defects that are due to faulty material or workmanship for a 2 year period from the original date of consumer purchase. If the product should become defective within the warranty period, we will repair any such defect without charge provided there is proof of purchase. This does not mean an extension of warranty period.

The warranty shall apply in all countries in which our product was sold via an authorized dealer.

This warranty does not cover the damage is due to by accident or misuse, abuse, alteration to the product or use inconsistent with the technical and /or safety instructions required.

This warranty shall not apply if the product has been dismantled or repaired by a person not authorized by us.

12. Copyright Reserved

Copyright © 2017 Qlarité SAS. All Rights Reserved.

Qlarité SAS reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Qlarité SAS. Qlarité SAS has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this documents. System design and specifications may change without prior notification.

Contact information for customer service:

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